otto schachnerPrinsessens Kvarter 2
7000 Fredericia, DenmarkEU Declaration of conformity available for download at
<http://www.os-safetycenter.com>

Batch no. 607-04-XXXXX.docx

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GB • User instructions**Brand name / type no.****Nilex Ultra PF**

607 070 Size 7
607 080 Size 8
607 090 Size 9
607 100 Size 10
607 110 Size 11

Description

Thin nitrile disposable glove with beaded cuff. The glove is powder free and ambidextrous.

Quality:

Weight: 4.0 g
Thickness of palm: 0.09 mm
Thickness of fingertips: 0.12 mm
AQL 1.5

Generally

Before use, it should be tested/ensured that the glove has the appropriate size to achieve the best possible comfort and safety at work.

The service life cannot be determined and depends on the scope of application and the extent to which the user makes sure that the glove is suitable for the intended use.

Category

The glove is certified in category III in compliance with Medical Device Regulation 2017/745 and PPE Regulation (EU) 2016/425 on safety requirements for personal protective equipment. The gloves were tested in accordance with the regulation no. 10/2011/EU and EU regulation no. 1935/2004 on Contact with foodstuff – Please check the specifications at the end of this document.

The glove has been tested in accordance with the standards EN 420:2003+A1:2009 (general requirements), EN 455 (medical use) and EN ISO 374-1:2016, EN 374-4:2013, EN ISO 374-5:2016 (chemicals and micro organisms).
Notified body responsible for EU Type Examination and Module C2 ongoing conformity assessment: af SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland (Notified Body No. 2777).

EN ISO 374-1:2016/Type B

KPT

EN ISO 374-5:2016

VIRUS

EN 374 Test Data / Type B**EN ISO 374-1:2016**

	Chemical	Time*	Level	Degradation %
A	Methanol	N/A	N/A	N/A
B	Acetone	N/A	N/A	N/A
C	Acetonitrile	N/A	N/A	N/A
D	Dichloromethane	N/A	N/A	N/A
E	Carbon disulphide	N/A	N/A	N/A
F	Toluene	N/A	N/A	N/A
G	Diethylamine	N/A	N/A	N/A
H	Tetrahydrofurane	N/A	N/A	N/A
I	Ethyl acetate	N/A	N/A	N/A
J	N-heptane	N/A	N/A	N/A
K	Sodium hydroxide 40%	>480	6	-42.9
L	Sulphuric acid 96%	N/A	N/A	N/A
M	Nitric acid 65 %	<10	0	97.6
N	Acetic acid 99 %	<10	0	93.9
O	Ammonium hydroxide 25%	>10	1	-52.0
P	Hydrogen peroxide 30%	>30	2	22.8
S	Hydrofluoric acid 40 %	N/A	N/A	N/A
T	Formaldehyde 37 %	>60	3	5.0
	**4% Chlorhexidine Digluconate	>480	6	19.0
	10-13% Sodium Hypochlorite	>480	6	14.7
	50% Sulphuric acid	>480	6	-20.5

EN 374-4:2013

10% Acetic acid	>120	4	66.7
5% Ethidium Bromide	>480	6	3.4
50% Glutaraldehyde	>480	6	27.4
0.1% Phenol	>480	6	33.8
1.5% Methanol in water	>480	6	21.9
70% Isopropanol	<10	0	62.2
35% Ethanol	<10	0	38.8
3% Povidone-iodine	>480	6	33.7
10% Sodium Percarbonate	>480	6	15.4

*) Breakthrough time (minutes)

**) Permeation rate: 7 µg/cm²/min

Performance Level	1	2	3	4	5	6
Minimum breakthrough time (mins)*	10	30	60	120	240	480

*) Glove performance quoted is based on laboratory data and may not reflect the actual duration of protection in the workplace due to other factors influencing the performance such as temperature, abrasion, degradation etc.

EN374-4:2013

Degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemical.

EN ISO 374-5:2016

Resistance to Bacteria and Fungi = Pass
Resistance to Virus = Pass

Migration test
(EN 1186-1:2002)**Application / qualities**

Thin and flexible seamless glove with great dexterity. The glove is waterproof and suitable for handling of many different chemicals with short contact time. If the glove has been in contact with substances it must be discarded, when the breakthrough time is reached. Primarily used in food production, mounting, product control, laboratories, chemical industry and for general cleaning.

Treatment / storage

The glove is for single use and should not be cleaned or re-used. Best stored in dark and cool surroundings in the original packaging.

Warning

This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals.
 The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only and relates only to the chemical tested. It can be different if the chemical is used in a mixture.

otto schachnerPrinsessens Kvarter 2
7000 Fredericia, DenmarkEU Declaration of conformity available for download at
<http://www.os-safetycenter.com>

Batch no. 607-04-XXXX.docx

- It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.
- When used, protective glove may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.
- Before usage, inspect the gloves for any defect or imperfections.
- The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimens.

Packaging100 gloves per dispenser.
10 dispensers in a carton made of recyclable cardboard.**SE • Bruksanvisning****Varumärke / art. nr.****Nilex Ultra PF**

607 070	Storlek 7
607 080	Storlek 8
607 090	Storlek 9
607 100	Storlek 10
607 110	Storlek 11

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Testdata re. Contact with Foodstuff: Regulation No. 10/2011/EU, EN1186, etc.:

overall migration for glove sample

Extraction Conditions	Analysis Result ; mg/dm ²	Permissible Limit ; mg/dm ²	Comment
3% (v/v) acetic acid aqueous solution @ 40°C for 1 hour	8.1	10	COMPLIED
10% (v/v) Ethanol aqueous solution @ 40°C for 1 hour	ND(<1)	10	COMPLIED
20% (v/v) Ethanol aqueous solution @ 40°C for 1 hour	ND(<1)	10	COMPLIED
50% (v/v) Ethanol aqueous solution @ 40°C for 1 hour	1.6	10	COMPLIED
Vegetable oil @ 40°C for 1 hour	ND(<2)	10	COMPLIED